



PRODUCT DATA SHEET

Product Code: ANI-W

High Level Disinfectant Wipes

- Powerful Disinfectant
- Tubs or Pouches
- Nonwoven Wipes
- Alcohol Free
- Effective Against a Range of Micro-Organisms
- Non-Corrosive

High level quaternary ammonium based disinfectant wipes. Manufactured from a nonwoven material, offering strength, absorbency and low lint properties at a cost effective price. Wipes are impregnated with a multi-active disinfectant solution that is alcohol free.

Recommended for the cleansing and disinfecting of hard surfaces of non-invasive, non-lumened medical devices such as monitors, IV poles and X-Ray Equipment.

These wipes are developed and tested by qualified chemists and microbiologists. They are then subjected to stringent EN testing procedures to ensure efficacy against a wide range of bacteria, viruses, fungi and spores, and are intended for cleaning and disinfection of hard surfaces.



V1.0 02.06.2021



PRODUCT DATA SHEET

Product Code: CHEM-W

This product has been tested effective against:

Virucidal Activity

- Adenovirus
- Hepatitis B
- Herpes Simplex
- Murine Norovirus
- Vaccinia Virus
- Enveloped viruses (including SARS CoV-2)

Bactericidal Activity

- Candida albicans
- Enterococcus hirae
- Escherichia coli
- Listeria monocytogenes
- Methicillin-resistant Staphylococcus aureus
- Methicillin-resistant Staphylococcus pseudintermedius
- Pseudomonas aeruginosa
- Salmonella tryphimurium
- Staphylococcus aureus

Compliance

This product has been tested against the following standards:

- EN 13727
- EN 13624
- EN 13697
- EN 16615
- EN 14476

V1.0 02.06.2021